**THE FREE GLIDING SCFE SCREW SYSTEM™**

**ENGLISH**

#### Device Description

The Free-Gliding SCFE Screw System is a self-extending cannulated screw intended for use in fixation of slipped capital femoral epiphysis and femoral neck fractures. The implant assembly includes a Male component (which is attached to the lateral cortex), a Female component (which anchors the femoral head) and a Cap. Anchorage of the components is achieved through threaded fixation. The screw has a built-in feature that allows for free extension of its length as the slipped capital physal plate heals and normal patient growth continues. Stable fixation and rotational stability is created at the fracture (slip) site while avoiding compression forces. The Free-Gliding SCFE Screw System is manufactured in medical grade Stainless Steel (316L, ASTM F138). It is available in 6.5mm and 7.3mm diameters. Two males (short and long) and 11 female sizes are combined to cover a range from 48 to 102mm in 2mm increments with a maximum of 30mm of growth potential.

#### Intended Use

The Free-Gliding SCFE Screw System is indicated as a temporary implant for stabilization of pediatric femoral neck fractures and slipped capital femoral epiphysis (SCFE) in all pediatric patients (less than or equal to 21 years old) with the exclusion of newborn and infants under 2 years of age.

#### Contraindications

Devices should not be used in patients with:

- Active or suspected latent infection or marked local inflammation in or about the affected area.
- Osteoporosis, insufficient quality or quantity of bone/soft tissue
- Compromised vascularity inhibiting adequate blood supply to the operative site.
- Material sensitivity documented or suspected.
- Sepsis
- Patients with abnormal neurological or mental conditions
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

Surgeons should warn patients about these contraindications and limitations when appropriate.

#### Adverse Effects

- Pain, discomfort, stiffness or abnormal sensations due to the presence of the device.
- Postoperative bone fracture and pain.
- Unrecognized joint penetration, chondrolysis
- Metal sensibility or allergic reaction to a foreign body.
- Bone resorption due to stress shielding.
- Premature physal closure, limb shortening, residual deformity
- Infection, both deep and superficial
- Avascular necrosis
- Nerve damage due to the surgical trauma
- Inadequate healing, continued slippage

#### Warnings

- Implants are single use items. Please note that single use device (SUD) which comes into contact with human blood or tissue should not be re-used and should be returned to the manufacturer or properly disposed.
- Metal implants should never be reimplanted. Although appearing undamaged, the device may have small defects or internal stresses which may eventually lead to implant failure.
- Correct implant handling is extremely important. Avoid contouring of metallic implants. Implant contouring or bending may reduce its fatigue strength and cause failure under load. Discard all damaged or mishandled implants.
- Failure to use largest possible components or improper positioning/insertion of the device during implantation can increase the possibility of loosening, bending, cracking, or fracture of the device or bone, or both.
- Continuous screening with an image intensifier (fluoroscopy) during guide wire insertion and whenever cannulated instruments are advanced over a guide wire is recommended to prevent unintended guide wire advancement and penetration into the surrounding tissues.
- Device breakage or damage can occur when implants are subjected to increased loading associated with delayed union, non-union, or incomplete healing.
- Implant contouring or bending may reduce its fatigue strength and cause failure under load.
- OrthoPediatrics Canada advises against the use of another manufacturer’s component with any OrthoPediatrics Canada component. Any such use will negate the responsibility of OrthoPediatrics Canada for the performance of the resulting mix.

**CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.**

#### MRI Safety Information

The Free-Gliding SCFE Screw System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of The Free-Gliding SCFE Screw System in the MR environment is unknown.

#### Surgical Technique

Standard surgical technique for the placement of the screws is recommended. Pre-operative procedures, knowledge of applicable surgical techniques, proper patient selection and correct implant placement of the implants are all equally important for the successful use of this product. The surgical technique manual details every step and should be carefully followed.

#### Device Lifetime and Retrieval

Removal or replacement of the implant is recommended subsequent to normal follow-up after closure of the physis or when full extension of the SCFE screw is observed. Routine removal of internal fixation devices may reduce the occurrence of symptomatic complications of implant breakage, implant loosening and implant related pain. In addition, if removal is favorable, OrthoPediatrics Canada recommends the retrieval of implants in order to avoid bone reduction and weakening, particularly in young and active patients. Although the final decision to retrieve the implants falls on the surgeon, a maximum Device Lifetime of 5 years for the implant has been defined to ensure material stability. The Surgical Technique manual details retrieval steps and should be carefully followed.

**Cleaning and Sterilization Instructions For Implant Components**
All implants are provided clean, but not sterile when shipped from OrthoPediatrics Canada. The instructions below should be followed for sterilizing items supplied non-sterile. Apply the standard cleaning protocol that is approved by the hospital before implant sterilization. All metallic implants can be steam sterilized following the instructions and parameters listed below.

**Note For USA only: Sterilization wraps, pouches, indicators and sterilization trays should be FDA-cleared for the sterilization cycle parameters.**

- Implant components of the Free-Gliding SCFE Screw system can be sterilized using sterilization pouches
- Devices should be dry before packaged for sterilization.

Method	Steam
Sterilization type	Prevacuum
Minimal temperature	270°F (132°C)
Minimal cycle time	4 minutes
Minimal drying time	60 minutes

#### Warning: Do not stack trays during sterilization

Other sterilization methods and cycles may also be suitable. However, validation of any alternative method using appropriate laboratory techniques is advised.

#### Cleaning, Sterilization and Re-sterilization Instructions for Instruments

Reusable instruments must be cleaned and sterilized prior to every use. The instrument tray and instruments of the Free-Gliding SCFE Screw System should be sterilized wrapped in two layers of 1-ply polypropylene wrap using sequential wrapping techniques. Please refer to document entitled: “Guidance for Instrument Care” for further information and instructions regarding cleaning, sterilization and re-sterilization of instruments.

#### Notice to the User and/or Patient

If you experience or are aware of any serious incidents that have occurred in relation to the device, please report them to OrthoPediatrics Canada at feedback@pegamedical.com. Additionally, kindly notify the competent authority of the Member State in which the user and/or patient is established. For Switzerland, please report them to materiovigilance@swissmedic.ch.

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